



# 6 in One

## SARS-CoV-2 & MP & RSV & ADV & Flu A/B

### Antigen Kit (Colloidal Gold)




## Features

- Reliable test performance
- Rapid result answers in only 15 minutes
- Easily visible and readable result
- Simple operation procedure, no equipment required
- Suitable for large-scale rapid screening
- 6 in One combo test. With only one swab delivers 6 results
- Quick differentiation of SARS-CoV-2, Influenza A/B, RSV, MP and ADV infection.

This Rapid Test kit is designed for health professionals and physicians in qualitative detection SARS-COV-2 / Influenza A / Influenza B / RSV / Adenovirus / Mycoplasma Pneumoniae in combo.

## Material provided

Sampling Swab



Test cartridge




Extraction Tube with buffer

# Specifications


Cat.No.	Pack Size
CG128001	1 Test/Kit
CG128025	25 Tests/Kit

## Operation




1

Collect sample



2

Mix liquid




3


Add sample then read results

## Result interpretation


### Possible Positive results



### Possible Negative results



### Invalid results



**Note:**

The test lines (A, B or T line region) may vary in intensity depending on the amount of antigen detected. Any visible (even light or faint) T line(s) should be interpreted as positive.

Co-infection with more than one pathogen is rare. If the results are positive for more than one antigen, e.g., influenza A, B or SARS-CoV-2, it's recommended to retest with a new sample.

Please refer the IFU for more details in section Interpretation of Results

## CONTACT

**Goldsite Diagnostics Inc.**  
 GOLDSITE Mansion,  
 No. 9 Rongtian Road, Kengzi,  
 Pingshan District, Shenzhen,  
 China, 518122

0086 755 2689 0807  
 export@goldsite.com.cn  
 en.goldsite.com.cn



@ Goldsite Diagnostics