


# SARS-CoV-2 & Influenza A/B Antigen Kit (Colloidal Gold)



## COVID-19 and Influenza

are both considered as highly contagious respiratory disease but triggered by different virus: SARS-CoV-2 and Flu A&B, respectively. GOLDSITE SARS-CoV-2 & Influenza A/B Antigen Kit (Colloidal Gold) helps healthcare professionals at the point of care to quickly diagnose and differentiate infections with any of the three respiratory viruses and aid in pursuing appropriate courses of action.




## Features


- Reliable test performance
- Rapid result answers in only 15 minutes
- Easily visible and readable result
- Simple operation procedure, no equipment required
- Suitable for large-scale rapid screening
- One swab delivers 3 results, quick differentiation of SARS-CoV-2 and Influenza A/B infection

## Material provided

Sampling Swab



Test cartridge




Extraction Tube with buffer

# Specifications

Cat.No.	Pack Size
CG125001	1 Test/Kit
CG125005	5 Tests/Kit
CG125025	25 Tests/Kit

# Operation



# Result interpretation

Results	SARS-CoV-2 results (Left window)	Influenza A/B results (Right window)
Positive		 Influenza A Positive
		 Influenza B Positive
		 Influenza A & B Positive
Negative		 Influenza Negative
Invalid		 Invalid Influenza test

## CONTACT

**Goldsite Diagnostics Inc.**  
 No. 103C, 503C & 504D, Technology Building & # 3A & 4A, Technology Building Annex, Zhaoshang Sub-District, Nanshan District, Shenzhen, China, 518067

0086 755 2689 0807  
 sales@goldsite.com.cn  
 en.goldsite.com.cn



@ Goldsite Diagnostics